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(GB) Operating Instruction
page 12-21

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(E) Manual de uso
página 32-41

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Liebe salivent® Kundin! Lieber salivent® Kunde!

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
Dear salivent® Customer!

Now you can easily utilise the healing power of the sea for your health and your wellbeing everywhere. salivent® assures effective inhalation for therapy and prophylaxis in the upper, middle and lower respiratory tract. Please read the following instructions thoroughly to achieve the best success with your salivent® Aerosol Inhalator: Breathe easily like being at the seaside.

Range of Applications:

The salivent® Aerosol Inhalator is used to inhale saline micro-fine aerosols in order to ease the breathing of the patients affected by the following illnesses:

- colds, such as head colds, coughs, hoarseness and the associated inflammation of the nose and throat
- asthma, bronchitis and other chronic illnesses of the upper and lower respiratory tracts
- allergies to house dust and pollen, as well as extreme dirt- and dust exposure
- heavy mucus- and crust formation in the nose and the associated obstruction of the nasal passages.

To achieve optimal results calmly breathe in and out during the whole inhalation process (inhale with the mask, exhale into the ambient air). This is important to enable the sole mist, i.e. the inhaled breathing air to reach the deep alveolar capillaries of the bronchia.

Loosened secretion should be coughed out by all means. You are free to take small breaks between individual inhalations and/or breaths.

6 Frequent Questions. 6 Answers

Why does the salivent® Aerosol Inhalator function so well?

It has been known for a long time that the inhalation of saline vapour (aerosol) supports the healing process actively. Aerosol particles are saline micro-particles, as they occur at the seaside at wind force 10. The inhaled aerosol stimulates the forming of mucus and secretion, simplifying coughing up and accelerating the natural self-cleaning mechanisms. A further positive effect of salt is its anti-inflammatory effect.

What are the holes and the rotatable rings for?

The holes control the air supply. The more holes are open, the easier it is to inhale. Convenient for persons with weak lungs such as children, the elderly or ill persons.

How long and how often should I inhale in order to feel better?

Depending on the illness we recommend to inhale at least 2-3 times per day. You may of course use the inhalator as often as you wish - overdose is not possible.

Is the salivent® Aerosol Inhalator allowed for use on children?

Yes, the device is suitable for children, as no noise is emitted, it is handy and can be used in any place. The mask expands to fit faces aged 4 and above.

Is it possible to inhale other salts?

Generally yes, as long as they are food-grade. For optimal healing however we recommend the salt provided with the inhalator. Grain size and type of salt are crucial conditions for 100 % effect.

Why can't I see any steam?

The inhalator is no electric evaporator from which steam evaporates but a self-active inhalation, i.e. the effect develops through the own power of the lung activity.

The Application

Step-by-Step Application!

Thoroughly clean your salivent® Aerosol Inhalator before every use also before the first use.

Fig. 1 Unscrew the screw cap from the housing.

Fig. 2 Now evenly sprinkle the grooves of the moistened (not dripping) sponge with the attached salivent® portion of sea salt (approx. 3g) to assure that the salt adheres well and is able to dissolve later.

Fig. 3 Alternative Dosage:
Fill the large container with water using the dosage aid included and the smaller container for the salt in the aluminum beaker, and then put the dry sponge onto the suction tube.

Fig. 4 Secure screw cap with the attached sponge onto the housing and close it hand tight.

Fig. 5 Attach the respiratory mask to the screw cap.

Fig. 6 Connect the salivent® Aerosol Inhalator to the mains supply, after 20 to 30 minutes the device warms up to ca. 38° Celsius and enables an agreeable inhalation at body temperature.

Fig. 7 Adjust the breathing holes by turning the valve ring to suit your personal needs and requirements.

Fig. 8 Now inhale for 30-45 minutes the soothing saline breathing air through your mouth and/or nose - both is possible. Breathe out into the ambient air.


1

2


3

4

5


6

20-30min

7

8

30-45 min

The Application

Important! General Information on Hygiene and Cleaning

We recommend the cleaning of the device after every use. For this, put the sponge, respiratory mask and screw cap for 1 to 2 minutes into boiling water. Bacterial germs and contamination will be killed in this way.

Ensure the individual-related use of the respiratory masks!

Fill up the housing with the aluminum beaker with boiling water and leave to set for 1-2 minutes before pouring it out. Do not rinse it under running water or placing it in the dishwasher. After use or cleaning of the device leave all components (especially the sponge) to dry.

NEVER PACK IT WHEN STILL WET!

Note!

During inhalation, salt crystals are formed in the respiratory mask, screw cap and aluminum beaker for physical reasons. This presents no risk to the user and no defect at the inhalator. These deposits can be easily removed during cleaning.

To ensure formation and thereby the effectiveness of the aerosols, we recommend the replacement of the application parts, in particular the sponges after approx. 100 applications.

If required you can order further respiratory masks and sponges from us.

We can also supply the salivent[®] sea salt which is provided in units of 20 portions.

Dosage and Duration of the Application

You can inhale as often and for as long as you wish, without any known disorders occurring.

However, to maintain and ensure the degree of efficiency we recommend:

1. to use a portion of salivent® sea salt (ca. 3 gr.) for every new application
2. to keep the sponge damp, not dripping wet, when preparing the application
3. to use the warmed-up device for approx. 30 minutes, (maximum 60 minutes)
4. to carry out an intermittent application, i.e. have breaks between inhalations
5. to carry out inhalations 2-3 daily over a period of 2-3 days

3 Things to Know

What else you should know

Due to the legal regulations for medical products we must point out the following:

1. The salivent[®] Aerosol Inhalator is only to be used for the purpose defined in the operating instructions.
2. For medical-technical equipment with electrical connection please note that in case of complaints the salivent[®] Aerosol Inhalator may only be repaired by the manufacturer or an authorised body.
3. We assume no liability for damages caused by improper use or non-observance of these instructions and/or by the use of parts and accessories which do not originate from our product range.
If you have a reason for complaint, please contact our customer service.

Warning

The salivent® Aerosol Inhalator complies with all relevant standards and guidelines in force at the time of production and selling. Our products have left the factory in a perfect condition. In order to maintain this condition the contents of this packing insert and the model plate on the inhalator must be complied with when handling this inhalator (transport, storage, installation, initial operation, application, servicing and closing-down), otherwise persons could be injured and the product and other equipment damaged. This safety-related information is only valid in the Federal Republic of Germany.

When used in other countries the relevant national rules are applicable. If the information in this packing insert is too limited, you can contact the manufacturer at all times under the address given.

- Please check the packing and complain immediately about any damage.
- Check before using the inhalator for the first time that it is not damaged. If in doubt, go to an electrical specialist for advice or contact the address given below.
- Before using for the first time check if the system voltage of the inhalator is identical with the mains supply of your country.
- Regularly check the inhalator for damage to the cable, plug, housing and other parts!
- If the inhalator is damaged, it is not to be used until it has been completely repaired!
- Repairs are only allowed to be carried out by authorised personnel. Please contact the address given below.
- Unplug it before cleaning! For external cleaning use a mild cleaning agent with a damp cloth. Never clean dripping wet!

3 Things to know

Technical Data:

Electrical Supply:	
Mains supply:	230 V AC
Mains frequency:	50 Hz
Power input:	max. 15 watt
Protection category:	IP 41
Protection class:	II

Performance Data:

Saline damp aerosol media throughput:	ca. 2.2 ml/min
Spectrum:	5-40 μm
Storage quantity:	200 ml
(measured with 0.9% NaCl in H ₂ O at 24°C)	

Aerosol Temperature

Temperature range:	37°C (+/- 3°C)
--------------------	----------------

Warning Signs

Signs and symbols (type plate on base of inhalator)

Read the operating instructions

Medical application part


No disposing of the inhalator with domestic waste

Device of Protection Class II. Protection isolated device.

Content

- 1 salivent® Aerosol Inhalator
- 2 respiratory masks
- 2 sponges
- 1 dosing aid
- 10 bags of salt
- 1 operation instruction

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Schäfer Produkte GmbH